


ADFS656



— ADFS656-0 - 0mm  
— ADFS656-100 - 100mm  
— ADFS656-50 - 50mm

波長 (nm)	0mm ...	50mm ...	100mm ...
300.00	0.217	0.012	0.003
310.00	0.330	0.019	0.006
320.00	0.258	0.014	0.005
330.00	0.113	0.006	0.002
340.00	0.071	0.003	0.001
350.00	0.122	0.007	0.002
360.00	0.662	0.036	0.010
370.00	5.374	0.272	0.081
380.00	22.644	1.076	0.323
390.00	43.044	1.965	0.590
400.00	51.570	2.327	0.697
410.00	53.047	2.397	0.718
420.00	53.294	2.417	0.724
430.00	53.492	2.429	0.732
440.00	53.598	2.436	0.734
450.00	53.672	2.448	0.737
460.00	53.757	2.459	0.743
470.00	53.841	2.462	0.741
480.00	53.970	2.474	0.747
490.00	54.091	2.483	0.752
500.00	54.231	2.486	0.752
510.00	54.332	2.501	0.755
520.00	54.462	2.513	0.759
530.00	54.552	2.525	0.763
540.00	54.677	2.526	0.764
550.00	54.755	2.539	0.768
560.00	54.913	2.545	0.772
570.00	55.015	2.562	0.777
580.00	55.157	2.568	0.777
590.00	55.164	2.571	0.782
600.00	55.340	2.582	0.788
610.00	55.424	2.586	0.789
620.00	55.604	2.598	0.787
630.00	55.779	2.606	0.793
640.00	55.932	2.631	0.801
650.00	56.010	2.636	0.803
660.00	56.030	2.645	0.803
670.00	56.193	2.638	0.808
680.00	56.308	2.648	0.804
690.00	56.303	2.668	0.816
700.00	56.487	2.672	0.812
710.00	56.723	2.691	0.826
720.00	56.969	2.724	0.832
730.00	57.355	2.762	0.826
740.00	57.450	2.778	0.832
750.00	57.606	2.804	0.842
760.00	57.752	2.814	0.851
770.00	57.759	2.827	0.853
780.00	57.762	2.837	0.859
790.00	57.855	2.860	0.868
800.00	57.883	2.890	0.874

[測定プロパティ]

波長範囲(nm): 300.00 ~ 800.00  
 スキャン速度: 高速  
 サンプルピッチ: 0.5  
 オート サンプルピッチ: 有効  
 測定モード: オート

[装置プロパティ]

装置の種類: UV-3600Plus シリーズ  
 測光値: 透過率  
 スリット幅: (20)  
 時定数: 0.1 秒  
 光源: 自動  
 光源切替波長: 290.00 nm

.....